

Optimized cell therapy manufacture process with Hydrogel-based CellFiber[®] encapsulation technology

Kevin Tran¹, Sindy Yang¹, Meriem Bejaoui², Selase Deletsu², Guang Han Ong², Nozomi Izumi², Tomoyuki Kaneiwa², Tianlun Lu², Iori Maeba², Li-Jie (Jerry) Wang²
 1. Business Development, CellFiber Co., Ltd., Tokyo, Japan. 2. Process Development & Innovation, CellFiber Co., Ltd., Tokyo, Japan.
 *Contact: kevin.tran@cellfiber.jp

Introduction

Regenerative medicine offers a powerful modality for treating cancer and tissue injury, yet its diverse indications and cell types introduce complex manufacturing challenges for clinical dosages.

CellFiber is a technology which encapsulates cells inside a hollow alginate hydrogel fiber for 3D cell culture (Figure 1). This fiber structure reduces exposure to hydrodynamic shear, maintains a defined cell-dense microenvironment, while remaining highly permeable to nutrients, oxygen, and metabolites. Here we present the integrated CellFiber culture platform and summarize its adaptation to iPSC, MSC, and T cell expansion, including growth kinetics, viability, and phenotype retention across cell types.



Figure 1. CellFiber[®] cell encapsulation technology

Method

[iPSC] Cryopreserved induced pluripotent stem cells (iPSCs) were thawed and pre-cultured under adherent conditions, then encapsulated and expanded for 7 days in hydrogel fibers within a 2 L Cellbag[™] on a Xuri[™] Cell Expansion System W25 (Cytiva).

[MSC] Human mesenchymal stem cells (MSCs) were thawed and pre-cultured under adherent conditions, then encapsulated and expanded for 8 days in hydrogel fibers within a 10 L Cellbag on the Xuri W25.

[T cells] CD3⁺ T cells were MACS-selected from cryopreserved PBMCs and encapsulated directly, without prior activation, then cultured for 12 days within a 2 L Cellbag on the Xuri W25. In a separate batch, CD3⁺ T cells were activated and transduced with LV-GFP at MOI = 2 at the time of encapsulation, then cultured for 6 days.

After expansion or transduction, cells were harvested by dissolving the hydrogel fibers in EDTA. Recovery, viability, and cell type-specific marker expression were assessed by flow cytometry and qRT-PCR (Figure 2).




Figure 2. CellFiber integrated cell manufacturing platform

Results

[iPSC]




Figure 3. iPSC expansion and differentiation capacity in fiber culture. (A) A yield of 8.8×10^9 cells (56-fold expansion) was recovered after 7 days at 97.9% viability (not shown), with pluripotency markers retained. (B) Following subsequent trilineage differentiation to germ layers, marker expression was higher in fiber-cultured cells than in 2D control.

[MSC]




Figure 4. MSC expansion in fiber culture and modeled process economics. (A) A yield of 2.2×10^9 cells (530-fold expansion) was recovered after 8 days at 91.6% viability (not shown), with ISCT-defined surface marker expression retained. (B) A cost model based on a projected batch yield of 1.10×10^{10} cells indicates up to 70% cost reduction versus conventional 2D multilayer flasks culture.

[T cells]




Figure 5. T cell expansion and transduction in fiber culture. A final yield of 2.0×10^{10} cells (750-fold expansion) was recovered at 92.2% viability after 12 days. Fiber culture yielded a higher proportion of CD4⁺ central memory and CD8⁺ naïve cells than 2D control (A). In a separate transduced batch, encapsulated T cells demonstrates significant increase of transduction efficiency with low lentiviral vector (MOI=2) without transduction enhancer (B).

Conclusion

- CellFiber[®] integrates encapsulation and expansion in a single closed process, and supports a combined encapsulation–activation–transduction workflow for T cells — reducing COGs, process length, and efficiency.
- Results on stem cells differentiation capability and T cells functional phenotype shows the advantages of CellFiber potentially improving therapeutic efficacy.
- A direct thaw-to-fiber workflow for MSC and iPSC, eliminating adherent pre-culture, is in development with availability targeted for late 2026.

